PL

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EN

## USER MANUAL BREATH ALCOHOL TESTER HT8G011

Thanks for purchasing our AAT088, a hand-held Digital Alcohol Breath Analyzer. AAT088 gives you a digital reading of your breath alcohol content and your possible concentration of alcohol intoxication.

### How Does Alcohol Testing Work?

When alcohol enters the blood through the alimentary canal (mouth, stomach, and intestines), it is distributed through every part of the body, including the lungs. As the blood passes through the lungs a quantity of alcohol passes through the alveoli (tiny sacs of air in the lungs). This quantity of alcohol that passes through the alveoli is proportionate to the amount of alcohol present in the blood. This period will last for about 15-30min before stable.

By using a breath testing device, it is possible to measure the quantity of alcohol that passes through the alveoli and consequently calculate the person's Blood Alcohol Concentration (known from here on as BAC).

The ratio of breath alcohol to BAC is 2000:1 - in simple terms, every 2000 milliliters of breath will contain the same amount of alcohol as 1 milliliter of blood. [and some other countries is using 2100:1 or 2200:1 or 2300:1]

For this reason, user need to blow for at least 5 seconds so that a sufficient deep air lung sample is obtained. Failure to do so may result in an artificially reading. Our device detects human breath and gives a digital reading out of a person's approximate BAC.

### Important

Don't be confused by what describes on internet. The BAC level in your body will change over time. It's not just a simple formula to calculate your weight, age and gender. Your body physical condition (tired, fasting & liver detoxification, etc...), will affect the speed of alcohol enter into your bloodstream. Your "BAC" will then gradually drops to zero when alcohol is absorbed and completely eliminated.

Typically, this period lasts from a few minutes to 10 hours, depending on how much you drink, your physical condition, and primarily your liver's ability to detoxify. If your liver is well detoxified, your BAC will be much lower than others

To begin, rinse your mouth with fresh water before using the device. Abstain from smoking or drinking for at least 15 minutes. If you perform a test just after drinking, the reading will fluctuate from very high to low

reading (i.e. it will drop gradually). It is due to the alcohol residues inside your mouth that effects the test values. This will also happen with aromatic drinks, alcoholic mouthwashes or mouth sprays, medicines or drug, as well as by burping and vomiting.

Please read the full instructions carefully for proper use of the AAT088 unit. Never drink and drive - the only safe level is zero!

## Auto cleaning and Battery Power

If first time of use or do not use more than 3 months, better perform Auto-cleaning to clear any oxide that cover sensor's surface. When power off, press & hold "ⓘ" & "Ⓞ" key for 3 seconds. It will show the voltage value of the battery and then "Lr" flashes. The machine will start heating up the sensor. User can press "ⓘ" to stop the auto cleaning at any time. After auto cleaning, battery power will be low and need to recharge the battery before normal use.


AAT088 comes with a built-in 3.7V Li ion battery (300mAH) and please fully charge the tester before use.

**ATTENTION:** The reading is invalid if the battery power is low.

When there is not enough battery power, the "Lo" & "bAt" will flash

Please connect the provided USB cable and charge the machine. The „Charge indicator“ will turn into red color. The red indicator will turn off when the battery is fully charged.

Insert the provided Mouthpiece into the gas inlet located aside of the LCD screen.

## INSTRUCTION

### - How to power on and start breathalyzer test

- (1) Press the Key to turn on the power and enter into the alcohol mode. "Wait" will turn on and activate sensor pre-heating mode, [ PH ] will flash for few times and then start a 15 sec. countdown. When countdown to zero, it is ready to start alcohol test.  
 (2) „Blow” flashes with beep tone and it is ready to take sample. There will have a 9 sec. count down period and you can either do 2.1 Manual Test or 2.2 Active Test

2.1 Manual Test. (if can not continuously blow for 5 second) Press the , key once, „Blow” lit on. Please start to blow through the mouthpiece steadily and continuously for approximately 5 seconds.

2.2 Active Test. Please blow in and "Blow" lit on. While blowing, a long beep will be heard and please keep blowing through the mouthpiece steadily and continuously for approximately 5 seconds. When the beep tone stops, sampling finishes and the test result reading will be shown on the display immediately.

The display will show animation while blowing as followin


If user do not contiguously blow for 5 sec. „Err” will show on display. Please press retest again and repeat (1)&(2).

Press , key can return back to standby mode

### -How to check test record

Under the "RLC" mode, press key for about 1 seconds to enter into the memory recall mode. Press , key again to scroll the records. (Max.10 records). Press and hold key for 2 seconds to exit memory recall mode.

In the record mode, the test number and its reading will flash on the LED alternatively. For example, the 9th record is 0.58%BAC (= 0.058%BAC), then r□g and 0.58%, will flash alternatively. If memories full, the new reading will overwrite the previous one (FIFO). All the records will be erased if user takes off the rechargeable battery.

### -How to power off:

(1) Auto off, the machine will turn off automatically if no key being press for about 25 seconds.

(2) Manual, press and hold the key for 1 second to power off the machine when display show RLC

### -How to select Unit/ Alert setting:

**Note: This is for service center use only, please do not change it unless you know well about your country regulation.**

(1) When display show "RLC" press "i" & "u" key at the same time for 3 seconds and "Unit" will display on the 8-segment display, press the "u" key to scroll 4 different units.

(2) After select the desired unit, either hold the "i" for 2 seconds and exit the setting mode, or press the "i" once to go to the alert setting mode.

Press the "u" key to scroll 4 different alert level **b-2, b-3, b-5, b-8** as following.

Alert value being set		Alert level (shown on display)			
		b-2	b-3	b-5	b-8
Unit	%oBAC	0.2‰	0.3‰	0.5‰	0.8‰
	%BAC	.02% (0.02%)	.03% (0.03%)	.05% (0.05%)	.08% (0.08%)
	mg/l	0.1 mg/l	0.15 mg/l	0.25 mg/l	0.4 mg/l
	mg/100mL	20 mg/100mL	30 mg/100mL	50 mg/100mL	80 mg/100mL

After select the unit and alert level, you can either leave the machine rest for 3 seconds to exit the UniUAlert setting mode automatically or hold the "i" key for 2 seconds and exit this setting mode.

### -How to read alert signals

AATOSS alcohol tester is designed to alert user as following:

Low Level	For measurements between 0-0.08 %oBAC (0-0.04mg/l), 0 will be displayed and avoid false alarm due to environment or oral smell in humans' mouth.
Safety Level	For measurements between Low Level and Alert Level, the alcohol concentration level shows on the display.
Alert Level	For measurements equal to or above the alert concentration level, the analyzer beeps together "HI" flashes for 10 times

### Caution

- (1) To achieve an accurate reading, user should ensure there is no alcohol residue inside the chamber before starting a new test
- (2) After 4 consecutive tests, please allow the device to rest for 2 min.
- (3) When the product is unused, please store it in a dry environment.
- (4) When low power, the reading is invalid, please recharge the machine immediately.
- (5) If the difference between two successive reading is >0.1%, user should test again.
- (6) If a very high concentration being tested, sensor's surface will be surrounded by lots of ethanol particles, which need to wait for few minutes for evaporation. Please wait for 3 min. for another test.

### Specification

Display Range: 0/ 0.08-2.00 o/oBAC or 0/ 0.008-0.200 %BAG or 0/ 8-200mg/100ml or 0/ 0.04-1.00 mg/l

(%o & mg/l conversion factor is about 2.0) Power Input: 3.7V (rechargeable battery).

Charger: 5V, 2A or above

Display: Digital LED display (3 decimal points) Record: 10 records

Operating Temperature: 5-c-5o-c, for optimum operation Agency Approval: CE UKCA FCC

Accessories: 8 mouthpieces, recharge battery and charge cable.

### **Precautions**

Your blood alcohol concentration may change within minutes Do not use this device to determine whether you can drive or not.

You cannot use the test results in Court.

Your safest bet is not to mix drinking with driving at all. Avoid Binge drinking.

### **Liability for Proper Function or Damage**

Never unscrew the Alcohol Tester. Warranty is not valid if any screws had been loosened or warranty label had been broken; or if the instrument is used in a manner not conforming to its intended use.

Other factors (e.g. environment, smoke & oral smell and other mouth impurities) may affect the BAC level, please rinse your mouth with fresh water before using the device. Abstain from smoking or drinking for at least 15 minutes and use the readings for Reference Only.

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